

# Acupuncture LASER Med-100 N

*Acupuncturist • Auriculo-Therapist  
General Practitioner • Dentist • Veterinarian • Neural Therapist  
Natural Practitioner • Physical Therapist*

In 1988 **LASOTRONIC** introduced the world's first Pocket-Therapy-Laser and only 5 years later - another world's first: The smallest, battery operated microprocessor controlled Acupuncture Laser **Med-100N** - small, beautiful and reliable. With its 130g weight (only 4.5oz!) it is a pleasure to hold, to operate and to work with. In this latest model, 7 menus lead you through its unlimited possibilities in auricular therapy, body laser acupuncture and any other form of LLLT - Laser Therapy.

Your last set parameters are being saved for your next need. This unique Laser is available at an affordable price and will secure your today's investment into the future. Choose between visible red and invisible infra red - choose from 20 to 40 mW output power (higher levels are not usually applied in laser acupuncture).



### Standard Accessories:

*Transportation case containing:*

- Conical treatment head
- Acupuncture glass fiber
- leather case for your pocket
- 1,5V Alkaline batteries
- integrated equalizer (power meter and timer)
- CE laser shields
- instructions manual

**Housing:** Aluminum, anodized in blue

### Options:

- Adapter to mains
- other glass fibers bended or straight

Model	Med-100 N red	Med-100 N infra red
<b>Wavelength</b>	660 nm (visible)	830 nm (invisible)
<b>Output Power</b>	<20 mW	<40mW
<b>Laser Diode</b>	GaAlInP	GaAlAs
<b>Laser Class</b>	3b	3b
<b>Weight</b>	130g	130g
<b>Dimensions</b>	172 x 32 x 16 mm	172 x 32 x 16 mm
<b>Input Voltage</b>	3 V DC	3 V DC
<b>Batteries</b>	2xAAA (Alcaline)	2xAAA (Alcaline)
Med. Prod. Class: IIa	CE 0297	ISO 9001 ISO13485 93/42EWG

**LASOTRONIC** reserves the right to change specifications in the interest of medical and technical developments without notification.

### 7 Menus: Double-quarz stabilized Med-100 N:

- 1) 7 Nogier-Frequencies: A:292 B:584 C:1168 D:2336 E:4672 F:9344 G:146 Hz + CW (Continuous Wave)
- 2) 7 Bahr-Frequencies: B1:599,50 B2:1199,0 B3:2398,0 B4:4796,0 B5:9592,0 B6:149,875 B7:299,75 Hz + CW
- 3) 10 Personal Frequencies (1Hz - 10 KHz) + CW
- 4) Nogier- & Bahr-Frequencies with timer function (15, 30sec, 1,2,3,4 & 5 min with automatic stop)
- 5) 10 Personal frequencies with timer function
- 6) Internal calibration (shows power in milliwatts)
- 7) Lock & Unlock (electronic safety lock)

**Info: Ask your local LASOTRONIC dealer or visit our homepage for:**

- *Bibliography:* a list of over 300 publications including double-blind-studies, clinical studies etc.
- *Indications (What is Laser Therapy):* a list of indications based on leading publications

**LASOTRONIC**<sup>®</sup>  
MEDICAL THERAPY LASERS

LASOTRONIC Intl. Ltd.  
Soroskari u. 48  
H-1095 Budapest

[www.lasotronic.com](http://www.lasotronic.com)